



Hugs® and Kisses® Approved Cleaning Materials

TECHNICAL BULLETIN



DESCRIPTION

This Technical Bulletin provides important information about which cleaning materials have been approved/not approved by Securitas Healthcare for use on Hugs Infant Protection and Kisses Mother-Infant Matching Tags.

Clients are advised to review the list of approved and non-approved cleaning materials carefully since manufacturers often provide a range of similar sounding product variants with different active ingredients. This may result in one variant of cleaning material being approved while another variant from the same manufacturer possibly not being approved. For example, Clorox Commercial Solutions® Germicidal Wipes are approved for use on Hugs and Kisses Tags whereas Clorox Healthcare™ H2O2 Cleaner disinfect wipes are not approved.

CLEANING INSTRUCTIONS

Hugs and Kisses Tags are not sterilized at the factory. They must be cleaned before first use and before inserting into tag charging station after patient use. In general, customers should consult with their Infection Control representative for cleaners available to the hospital that have been approved for use on plastics and are on the Approved Cleaning Materials list in this Technical Bulletin.

Customers must follow all applicable label instructions for these products. Many of these products specify a two-step process with specific surface contact times and/or instructions for rinsing away residue remaining on the tag.

Note: If residue remains on the tag, the water evaporates away, leaving the chemicals to increase in concentration. A normally benign, mild alkaline solution then becomes a strong (and damaging) alkaline solution. This could cause serious degradation to a tag's plastic casing. Refer to the MSDS for cleaning specifications and dry times.

If required, a disposable soft-bristle brush may be used to remove surface dirt from the tag's surface and to gently clean inside the tag prongs.

For minimal surface degradation and cracking over time, use the following guidelines:

- Do not use pathogenic cleaners specified for tuberculosis (TB).
- Do not soak the tag.
- Do not use an autoclave to clean the tag or serious damage may result.



APPROVED CLEANING MATERIALS

The following list of cleaning materials has been approved by Securitas Healthcare for use on Hugs and Kisses Tags:

Cleaning Materials (Approved)

3M™ HB Quat Disinfectant Cleaner Concentrate 25
3M™ Neutral Quat Disinfectant Cleaner Concentrate
Brulin® Maxima 128
CaviCide™ or CaviWipes™
Clorox Commercial Solutions® Germicidal Wipes
Clorox® germicidal bleach (bleach use with gauze)
Coverage™ HBV
Ecolab® Airkem® A-456-N
Ecolab® Asepti-zyme™
Ecolab® Optipro™ Solid Multi Enzymatic Manual Cleaner
Expose® 256
Expose® II 256 or Phenolic 256 DC™
Formula C™
Hibiclens
Frekaderm Actief 100 g liquide
Lingettes Prodene Aseptil (FDSUL233)
Metrex CaviWipes™
Metrex Detergezime™ - generic detergent + enzyme
Oxivir® Five 16 Smart Dose
PDI® Super Sani-Cloth® germicidal disposable wipes (purple top)
Precise® Hospital Foam Cleaner Disinfectant
Premoistened Alcohol/DI Clean Wipes
Professional Amphyl® Hospital Bulk Disinfectant
Prolystica® Ultra Concentrate Enzymatic Cleaner
Sani 24
Sani-Cloth® AF3 Germicidal Disposable Wipes
Sani-Cloth Bleach Orange Top
SaniWipes / Red Top = Sani-Cloth® Plus

Spor Klenz® - Ready To Use
Staph-Attack™
Thymo-Cide
Vesphene® II SE
Virkon® -S Powder
Virox™ Accel® Prevention (wipes & concentrate in use)
Virox™ Rescue™
Wexford Labs Ready-to-use Wex-cide TM germicidal detergent
Wexford Labs Wex-Cide 128 concentrate

NON-APPROVED CLEANING MATERIALS

The following list of cleaning materials has not been approved by Securitas Healthcare for use on Hugs and Kisses Tags. **DO NOT USE** these cleaning materials:

Cleaning Materials (Non-Approved)

Azul® Staphicide -Eco Line 0.5%
Bulk Disinfectant - not defined - many versions
Clorox Healthcare™ H2O2 Cleaner disinfect wipes
FDS Surfianos Premium NPC
Hydrogen Peroxide Cassettes
Isopropyl Alcohol use with gauze
LpH® Disinfectant Cleaner
Metrizyme®
Oxivir® TB/Oxivir® TB wipes
Sporicidin®
Virex® II 256
Wexcide 128
Wexcide-Ready-To-Use

4600 Vine Street
Lincoln, NE 68503
+1 888 622 6992

North America
info@securitashealthcare.com

International
global@securitashealthcare.com

securitashealthcare.com

About Securitas Healthcare

Securitas Healthcare empowers caregivers to deliver connected, productive and safe care. Its innovative portfolio of healthcare solutions helps over 15,000 hospitals, clinics and senior living organizations worldwide protect people, use assets efficiently and understand their operations for a caring and healing environment. Securitas Healthcare is proud to be part of Securitas, the world's leading intelligent protective services partner. For more information, visit us at securitashealthcare.com

Doc: 0970-090-000 REV H | Published: 2024/01/10 | KB Article 8184

Hugs and Kisses are trademarks of Securitas Healthcare LLC and/or its affiliates.
© 2024 Securitas Healthcare LLC. All rights reserved. This publication is accurate as of the date of printing.