

M210 Fall Monitoring System

SOLUTION OVERVIEW



Supporting safety and individualized care

SOLUTION BENEFITS

- Address many needs with one solution. No need to purchase and maintain multiple solutions to meet the diverse needs of the individuals in your care. Choose the delay times, tamper resistance and alert type settings that best support each person's situation and care plan.
- Give caregivers more time to focus on people. The M210 reduces the need to continually disconnect and reconnect pad cords. It helps reduce nuisance alarms and alarm fatigue, and enables caregivers to spend more time with patients or residents.
- Simplify maintenance and operation. The M210 firmware can be updated with its mobile app, eliminating the need to regularly replace monitors. Standardizing on a single solution reduces training requirements and streamlines ongoing maintenance and operation.

CHALLENGES

Implementing and maintaining a fall management program to reduce the quantity and severity of falls. Developing and supporting individual care plans, including fall prevention measures.

SOLUTION


The M210 Fall Monitoring system is an enhanced version of our M200 fall monitors, trusted by healthcare facilities for over 20 years. Sleek and durable, the multi-featured M210 offers a single, easy-to-configure platform to support your unique fall management program.

HOW IT WORKS

With the M210 Fall Management System, you can connect two pressure-sensitive pads to a single monitor (one on a bed, another on a side chair or wheelchair). Select a time delay setting that best protects each person in your care. When a patient starts to exit the bed or chair, the monitor issues visual and auditory alarms to alert staff that immediate assistance is required. If a nurse call system is connected, alarms can be sent to the nurse call station.

The M210 monitor is easy to configure with a mobile app, allowing you to wirelessly configure the fall monitor to the individual needs of each patient or resident. Choose the configuration that best suits the criteria to successfully execute care plans and program goals: electric or battery power, auto-reset or non-auto-reset, audible or silent alerts, and tamper-resistant or easy-access rear cover.





Simplify configuration and deployment

The M210 Fall Monitor can be easily configured and its firmware updated using the mobile app.

4600 Vine Street Lincoln, NE 68503 +1888 622 6992

North America info@securitashealthcare.com

International global@securitashealthcare.com

Long-term partnerships, long-term value

Exceptional care requires teamwork, which is why we select only quality partners to deliver great products and service at the best price. Each of our trusted providers offer an unmatched level of personalized service and support to ensure end-to-end success. We value our partners and further stand by them with additional training to ensure excellent client service and support at each stage of implementation.

About Securitas Healthcare

Securitas Healthcare empowers caregivers to deliver connected, productive and safe care. Its innovative portfolio of healthcare solutions helps over 15,000 hospitals, clinics and senior living organizations worldwide protect people, use assets efficiently and understand their operations for a caring and healing environment. Securitas Healthcare is proud to be part of Securitas, the world's leading intelligent protective services partner. For more information, visit us at securitashealthcare.com

DOC-12-70036-AG